

Labtex

AATCC Tumble Dryer For Home Laundering Test

Model: **LBT-M6D**

Meet Standard

**AATCC 88C, I24, I35, I43,
I50, I72, I79, I88**

(AATCC Laboratory Procedure I-2018)

Cycle	Normal Permanent Press, Delicate
Capacity	15KG
Power	220V/50Hz OR 120V/60Hz
Packing Size	683mm x 711mm x 1092mm
Packing Weight	61KG



AATCC Standard Tumble Dryer Parameters

Cycle	Normal	Permanent Press	Delicate
Exhaust Temperature	68±6°C	68±6°C	60±6°C
Cool Down Time (Min)	≤10	≤10	≤10