

Labtex

AATCC Washing Machine For Home Laundering Test

Model: **LBT-M6**

Washing Cycle	1 For AATCC LP1 - Normal - Permanent Press, - Delicate 2 Drain & Spin 3 Rinse & Spin 4 Customized*
Water Control	25°C~60°C±1°C (wash) Tap water (Rinse)
Washing Load	10.5KG
Power	220V/50Hz/1KW OR 120V/60Hz/1.5KW
Packing Size	700mm x 680mm x 1040mm
Packing Weight	78KG

* Base on customer needs to tailor made the specific washing sequence with LBT-Portable Handheld Programmer.



OPTIONAL:

LBT-Portable
Handheld programmer



Meet Standard

**AATCC LP1, 88B, 88C, 124,
130, 135, 142, 143, 150, 159,
172, 179**

Washing Cycle for AATCC LP1-2018

Cycle (Table I)	Normal	Permanent Press	Delicate
Water Level Medium	19±1gal	19±1gal	19±1gal
Agitation Speed	86±2spm	86±2spm	27±2spm
Washing Time	16±1min	12±1min	8.5±1min
Final Spin Speed	660±15rpm	500±15rpm	500±15rpm
Final Spin Time	5±1min	5±1min	5±1min

Water Temperature Control

Wash	Rinse
From 25°C to 60°C±1°C. Setting Washing Temperature in 27°C, 41°C 49°C, 60°C for AATCC Standard, or any other temperature will be specified.	Tap Water