

DRK107A Paperboard Thickness Tester



Product Introduction

DRK107A Paperboard thickness tester is a special instrument for measuring corrugated board.

Product Features

Manual type, digital display/pointer type and percentile/micrometer type, compact and portable structure.

Product Application

This equipment is suitable for the thickness measurement of flat-sheet pattern and corrugated cardboard.

Technical Standard

Product parameters

Item	Parameter
measuring range	0-25.4mm
Contact Area	
Contact Pressure	
Scale Division Value	0.001mm
Repeatability	$\pm 2.5\mu\text{m}$ (or $\pm 0.5\%$)
Dimension	240×160×120(mm)
Weight	4.5 kg

Main fixtures

- ✓ Mainframe
- ✓ One operation manual.

