

DRK113-1 RCT Sample Cutter



Product introduction

The instrument according to GB/1048 regulations to manufacture, it is the sample cutter for the test of all kinds of paper and paperboard.

Product features

The sampling devices pleasing to look, compact structure, reasonable and effortless operation, easy to use

Product application

It is the sample cutter for the test of all kinds of paper and paperboard. It is the supporting tool of DRK113 Crush Tester.

Technical standards

The instrument according to GB/1048 regulations to manufacture.

Product parameter

1、 Sample width: 12.7-0.025mm

2、 Sample length: 152-0.25mm

3、 sample both sides parallelism: $<0.015\text{mm}$

Main fixtures

Mainframe; operating manual; certificate of quality