

## DRK113-2 ECT sample cutter

### product introduction

ECT sample cutter is specific sample cutter for corrugated cardboard for ECT、PAT, it is supporting tools of DRK113 crush tester.



### product feature

Structure is simple 、 high accuracy 、 easy maintenance

### product application

The instrument is used for cutting sample for all kinds of cardboard, available for 3 layer 5layer、 7layer etc

### technical standard

This instrument is designed and manufactured according to the specified requirements of GB6543 and GB6548

### product parameter

- 1、 Verticality of blade movement line and the baffle:  $\pm 0.10\text{mm}$
- 2、 Depth of parallelism of blade of movement line and the locating end face of the locating board:  $\leq 0.10\text{mm}$
- 3、 Size tolerance of the cut sample:  $\pm 0.5\text{mm}$

### product configuration

mainframe、 operational manual