

DRK113-3 FCT sample cutter

Product introduction

FCT sample cutter is the Special apparatus of my company produced corrugated board flat crush test, it is the supporting tool of DRK113 Crush Tester.

Product features

The sampling devices pleasing to look, compact structure, high accuracy easy maintains.



Product application

Take a corrugated board, press the sampler in the corrugated board, Press the handle clockwise rotation to 3-5 laps and sample will be cut. Do not counterclockwise rotation to prevent the blade being damaged.

Technical standards

ISO3035 《Corrugated fiberboard — Determination of flat crush resistance》

GB/T 2679.6 《Determination of corrugated paper flat crush resistance》

Product parameter

Sampling area: 64.5 cm² 32.2 cm²

Main fixtures

Mainframe; operating manual; certificate of quality.