

## DRK114C Circle Sample Cutter



### Introduction

DRK114C Circle sample cutter is the sample cutter for the test of all kinds of paper and paperboard. The cutter conforms to the standard of QB/T1671—98.

### Characteristics

The instrument simpler and smaller, it can quickly and accurately cut the standard area about 100 square centimeters.

### Application

DRK114C Circle sample cutter is dedicated sampling devices for paper and paperboard physical performance testing, it can quickly and accurately cut the standard area about 100cm<sup>2</sup>.

## Standards

The instrument conforms to the standards of GB/T451, ASTM D646, JIS P8124, QB / T 1671.

## Parameter

Items	Parameter
Specimen Area	100cm <sup>2</sup>
Specimen Area error	±0.35cm <sup>2</sup>
Specimen thickness	(0.1~1.5)mm
Dimension Size	(L×W×H)480×380×430mm

## Configurations

Mainframe; Certificate of quality